

# OUR COMMITMENT TO YOU, IN THESE TIMES OF CHANGE

Covid-19 Prevention Plan

## Information about the Covid-19 tests distributed at the Grande Real Santa Eulália Resort & Hotel Spa, Real Marina Hotel & Spa and Real Marina Residence, during June holidays:

These are rapid tests for the qualitative detection of the SARS-COV-2 antigen (research of virus surface proteins by chromatographic immunoassay).

Brand: Roche Laboratories

Instructions for use: Insert the swab about 2 cm in the nasal cavity in both nostrils for approximately 15 seconds. After collecting the biological material, place the swab in the extraction tube, squeeze at the bottom and turn it about 10x. Close the tube with the dispensing cap. Place the 4 drops on the test plate, in the round field marked.

Read the result in 15-30 min after placing the 4 drops (the Control line - zone C will always appear, if a 2nd line appears on the T, it means that the test is positive).

They are not as evasive as PCRs (molecular test of viral RNA research exclusively performed in laboratories) where the collection of material in the swab is done more deeply in the nasopharyngeal cavity. This is exactly why these rapid tests are safe self-test tests without the need for a prescription.

Note: Real Hotels Group will not be held responsible for the results of the tests or for any infections that may occur during stays.

